

# Body Mold Girdles

---



## Medium Control Panty Girdle

#2000

- Extra High Waist



## Above-Knee Girdle

#5000

- Extra High Waist
- Side Zippers with Hook-Eye Closures
- Detachable Crotch Piece



## Below-Knee Girdle

#7000

- Extra High Waist
- Side Zippers with Hook-Eye Closures
- Detachable Crotch Piece



## Ankle-Length Girdle

#9000

- Extra High Waist
- Side Zippers with Hook-Eye Closures
- Detachable Crotch Piece

---

Please refer to the garment sizing chart for appropriate size

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Body Mold Full-Body Girdles

---



**Full-body Panty Girdle**  
#2700



**Full-body Above-Knee Girdle**  
#5700



**Full-body Below-Knee Girdle**  
#7700

- Side Zippers with Hook-Eye Closures -
- Adjustable Shoulder Straps -
- Detachable Crotch Piece -



**Full-body Ankle-Length Girdle**  
#9700

Please refer to the garment sizing chart for appropriate size

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.  
7758 Hampton Place  
Loganville, GA 30052  
800.525.7752  
770.934.0235  
[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Bras, Vests and Stabilizers



## Surgical Bra

#1524

Available in sizes from 34"-50"

- All Stretch Cup
- Front Closure
- Adjustable at Shoulder



## Compression Bra

#700 Fits A/B - 30-32

#701 Fits C/D - 34-36

#702 Fits D+ - 38-40

- Same Compression Fabric as our Girdles
- Front Closure
- Adjustable at Shoulders



## Sub-Pectoral Implant Stabilizer

#300



## Female Vest

#4100

- Adjustable at Shoulders

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Bands and Binders

---



#1100  
9" Abdominal Band  
S - 22-26"  
M - 26-30"  
L - 30-34"  
XL - 34-38"  
XXL - 38-42"



#1150  
12" Abdominal Band  
S - 22-26"  
M - 26-30"  
L - 30-34"  
XL - 34-38"  
XXL - 38-42"



#1200  
3- Panel Binder  
30-45"

#1200XL  
3- Panel Binder  
46-60"



#1250  
4- Panel Binder  
30-45"

#1250XL  
4- Panel Binder  
46-60"

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)



# Body Mold Facial Support and Compression Sleeves

---



**Facial Chin/Neck Support**  
#1000



**Compression Sleeves**  
#4500

-Adjustable front and back

**Please refer to the garment sizing chart for appropriate size**

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Male Compression Garments



**Male Support Briefs**  
#2100



**Gynecomastia Vest**  
#3100



**Full Body Men's Compression Girdle**  
#3700

Purchase from:



Cosmetic Surgery Suppliers, Inc.  
7758 Hampton Place  
Loganville, GA 30052  
800.525.7752  
770.934.0235  
[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Pic-Eze™

---

- **Pic-Eze** are disposable fabric panties that can be worn for pre and post operative photos, exams and many aesthetic procedures.
- **Pic-Eze** are made from a soft latex-free material. Each individually packaged **Pic-Eze** provides patients comfort and dignity at a low cost to the surgeon.



## Pic-Eze™

- Disposable
- Sanitary
- Comfortable
- Economical
- Breathable Mesh
- Flexible
- Latex Free

- **Pic-Eze** come in packages of 10 individually packaged panties and are available in 2 sizes.

<u>Size</u>	<u>Waist</u>
Large	20 - 36
X-Large	26 - 48

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Garment Sizing Chart

## Female Girdles [Cat.# 5000, 7000, 9000, 2700, 5700, 7700, 9700]

	Waist	Hips
Small	<b>23-25"</b> (58.5-63.5cm)	<b>28-36"</b> (71.0-91.5cm)
Medium	<b>25-29"</b> (63.5-73.5cm)	<b>30-39"</b> (76.0-99.0cm)
Large	<b>29-33"</b> (79.5-84.0cm)	<b>34-43"</b> (86.5-109.0cm)
X-Large	<b>33-37"</b> (84.0-94.0cm)	<b>40-49"</b> (101.5-124.5cm)
XX-Large	<b>37-41"</b> (94.0-104.0cm)	<b>47-57"</b> (119.5-145.0cm)
XXX-Large	<b>41-46"</b> (98.4-110.4cm)	<b>55-65"</b> (132.0-156.0cm)

## Arm Sleeves [Cat.# 4500]

	Bicep
Small	<b>9<sup>1/4</sup>-10"</b> (23.5-25.5cm)
Medium	<b>10-11<sup>1/2</sup>"</b> (25.5-29.0cm)
Large	<b>11<sup>1/2</sup>-13"</b> (29.0-33.0cm)
X-Large	<b>13-15"</b> (31.5-36.0cm)

## Female Vest [Cat.# 4100]

	Bra
Small	<b>34-36"</b> (86.5-91.5cm)
Medium	<b>36-38"</b> (91.5-96.5cm)
Large	<b>40-42"</b> (101.5-106.0cm)
X-Large	<b>44-46"</b> (111.0-116.0cm)
XX-Large	<b>48-50"</b> (115.2-120 cm)

## Abdominal Band [Cat.# 1100-1150]

	Waist
Small	<b>22-26"</b> (56.0-66.0cm)
Medium	<b>26-30"</b> (66.0-76.0cm)
Large	<b>30-34"</b> (76.0-86.5cm)
X-Large	<b>34-38"</b> (86.4-96.5cm)
XX-Large	<b>38-42"</b> (91.2-100.8cm)

## Abdominal Binder [Cat.# 1200-1250]

	Waist
1200	<b>30-45"</b>
1250	<b>30-45"</b>
1200XL	<b>46-60"</b>
1250XL	<b>46-60"</b>

## Male Vest [Cat.# 3100]

	Chest
Small	<b>36-38"</b> (91.5-96.5cm)
Medium	<b>40-42"</b> (101.5-106.5cm)
Large	<b>44-46"</b> (112.0-116.5cm)
X-Large	<b>48-50"</b> (115.2-120cm)
XX-Large	<b>52-54"</b> (124.8-129.6cm)

## Male Girdle [Cat.# 2100]

	Waist
Small	<b>28-32"</b> (71.0-81.5cm)
Medium	<b>33-36"</b> (84.0-91.5cm)
Large	<b>37-40"</b> (94.0-101.5cm)
X-Large	<b>41-44"</b> (104.0-112.0cm)
XX-Large	<b>45-49"</b> (108-117.6cm)

## Full Body Male [Cat.# 3700]

	Waist	Chest
Small	<b>28-32"</b> (71.0-81.5cm)	<b>36-38"</b> (91.5-96.5cm)
Medium	<b>33-36"</b> (84.0-91.5cm)	<b>40-42"</b> (101.5-106.5cm)
Large	<b>37-40"</b> (94.0-101.5cm)	<b>44-46"</b> (112.0-117.0cm)
X-Large	<b>41-44"</b> (104.0-112.0cm)	<b>48-50"</b> (115.2-120cm)
XX-Large	<b>45-49"</b> (108.0-117.6cm)	<b>52-54"</b> (124.8-129.6cm)

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)



# Liposuction Aspirator

---



---

## Our liposuction aspirator features:

- Advanced state of the art piston pumps
- Maximum deep vacuum quickly reached and maintained
- Dual headed pumps eliminating back pressure
- Pumps set in series rather than in tandem
- Three-quarter horsepower
- Quiet, vibration-free, reliable operation
- Standard implosion-proof canisters and disposable liners
- Five-year warranty on all parts and labor

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

**cosmeticsurgerysuppliers.com**

# Aspirator Accessories

---



Standard Suction Tubing



Flexi-Tubing



Disposable HEPA Filter



2000cc Abbott Canister  
with Disposable Liner



1500cc Baxter Canister  
with Disposable Liner

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752



770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Syringe Aspiration

---


syringe cannula for toomey tip  
(made to order)



toomey tip syringe

syringe-to-syringe transfers

---



## Syringe Locks

---



snapper type

- 10 cc
- 30 cc
- 60 cc



horseshoe type

- 10 cc
- 30 cc
- 60 cc

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Disposable Tubing

---



Standard Tubing

- 3/8" I.D.
- Transparent for easy visualization
- 10 foot lengths for easy use
- 10 pieces to a case



Flexi-Tubing

- 9/32" I.D.
- Super flexible tubing
- 12 foot lengths for easy use
- 10 pieces to a case

## Tumescent Infiltration Tubing

---



AFT Tubing

- Bag spike and luer-lock connector
- Thicker section is silicone, to avoid latex reactions
- 10 pieces to a case



T.I. Tubing

- Bag spike and luer-lock connector
- 8 foot lengths for convenience
- 10 pieces to a case

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)



# Cannula & Brushes

**Cosmetic Surgery Suppliers, Inc.** Cannulas are custom made to your specifications when you place your order. When ordering your custom cannula please specify the following four items: tip style, lumen diameter, lumen length and handle style. All of the most popular tip styles are available and custom tip styles can also be manufactured.



One-Hole Standard



Three-Hole Mercedes



Three-Hole Accelerator



Candy Cane




Spatula




Three-Hole Cage

Standard Diameters  
2mm, 3mm, 4mm, 5mm & 6mm


Standard Lengths  
10cm, 15cm, 23cm & 30cm



One-piece cannula handle



Luer lock hub



60cc syringe hub

## Cannula Cleaning Brushes

Available in 5 sizes ( 2mm, 3mm, 4mm, 5mm & 6mm ), the brushes fit the full range of cannula sizes currently being used in most surgeries



**Cosmetic Surgery Suppliers Inc.** Cannula and tumescent needles are manufactured to the highest industry standards for quality. They include a two year warranty on materials & workmanship

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Tumescent Infiltration

## The Fluid Transfer Kit

### Kit contains:

- One sterile transfer system with dual-check valve and bag spike
- One sterile 60cc syringe
- Priced lower than current systems

For use where ever closed fluid transfer is desired. Large 60cc syringe allows rapid movement of fluids with fewer repetitions. Dual-check valve means no disconnect/reconnect problems.



## Blood Pressure Cuff

- Pressure infuser
- Fast fluid level inspection
- Durable and easy to clean
- Precise easily read pressure gauge



## Variable Speed Pump

- Flow rate controlled and adjustable from 35 - 350 ml/min
- Easy set up. Disposable tubing set
- Continuous infusion
- Pneumatic footswitch included
- Directional switch - tubing can be installed in either direction



Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgesuppliers.com](http://cosmeticsurgesuppliers.com)

# Tumescent Instruments and Tubing



Infusion Handle



Infusion Handle with On/Off Control



## Tumescent Needle

Luer Lock Connector With Multi-Directional Holes For Rapid Infusion of  
Tumescent Fluids

### Standard Sizes

2mm x 15cm

2mm x 23cm

2mm x 30cm

2.5mm x 15cm

2.5mm x 23cm

2.5mm x 30cm

3mm x 15cm

3mm x 23cm

3mm x 30cm



### AFT Tubing

(For use with electric pumps)



### T.I. Tubing

(For use with pressure bag)

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)



# Lipo-Injection Products

---



Reinjection Gun:

Available 1, 3, 5 & 10cc - 3, 5 & 10cc - 35cc - 60cc



UM-14 Underminer / Suture Puller




Type 1 Reinjector




Type 2 Reinjector



Type 3 Reinjector



Syringe-To-Syringe  
Transfer



Syringe Rack (autoclavable)  
Available 10cc, 35cc, 60cc

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)



# Basic Lipoinjection Kit

---



This kit contains the basic tools for lipo-harvesting and reinjection.

The kit is affordably priced and contains:

- 1 - 60cc harvesting cannula
- 1 - small tumescent needle
- 1 - 60cc syringe snap lock
- 1 - 14 gauge underminor with blunt needle
- 1 - 10cc reinjection gun
- 1 - 16gauge x 6cm type II lipoinjection needle
- 1 - 16gauge x 9cm type II lipoinjection needle
- 1 - 16gauge x 9cm type II angled lipoinjection needle

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.  
7758 Hampton Place  
Loganville, GA 30052  
800.525.7752  
770.934.0235

**cosmeticsurgerysuppliers.com**

# Fat Graft/Transfer in any Volume



## Tissu-Trans® FILTRON systems

**3-TT-FILTRON 500cc**


**3-TT-FILTRON 2000cc**

Simple, sterile non-traumatic harvest and filter of lipocytes.  
For all volume procedures:

- Eyelids
- Lips
- Face
- Hands
- Body Contouring
- Buttocks
- Breast Reconstruction
- Breast Augmentation
- Correction of Lipoplasty Deformities

Instructional Videos online:  
[www.shippertmedical.com](http://www.shippertmedical.com)

The Tissu-Trans® FILTRON systems do not change the normal flow of the physician's technique. Standard tumescent solution is infused. The Tissu-Trans® is placed between the harvest cannula and the waste canister. During harvest, the aspirate flows into the filter within the canister. The waste (blood, epinephrine, Xylocaine) is filtered and suctioned out, leaving quality fat that is ready for re-injection.



Open-ended, beveled STRAW cannula for removal/transfer from FILTRON systems.

## Fat Graft Instrument and Cannula Guide

Recommended equipment for autologous fat transfer/fat grafting by procedure.

PROCEDURE Re-injection Site	HARVEST Instruments		RE-INJECTION Instruments	
	Tissu-Trans® Model	Harvest Cannula with Biplane Handle	Re-injection Syringe	Re-injection Cannula
<b>EYELIDS/LIPS</b> up to 500cc	3-TT-FILTRON 500	3-H-SB-24H-3x30 24hole, 3mm x 30cm	3-S-LL1 1cc Luer Lock Syringes	3-C-1H2-18x8 1hole, 18gauge x 8cm
<b>FACE/BODY</b> up to 500cc	3-TT-FILTRON 500 face	3-H-SB-24H-3x30 24hole, 3mm x 30cm	3-S-LL5 5cc Luer Lock Syringes	3-C-1H2-16x8 1hole, 16gauge x 8cm
	3-TT-FILTRON 500 body	3-H-SB-4H-3x30 4hole, 3mm x 30cm	3-S-LL10 10cc Luer Lock Syringes	3-C-1H2-14x8 1 hole, 14gauge x 8cm
<b>BODY</b> Breast/Buttocks up to 2000cc	3-TT-FILTRON 2000	3-H-SB-24H-3x30 or 3-H-SB-4H-3x30	3-S-LL60 60cc Luer Lock Syringes	3-CTT-1H2-3x15 1hole, 3mm x 15cm
<b>BODY</b> Breast/Buttocks up to 1500cc	3-TT-MEGA 1500 Decantation Method	3-H-SB-24H-3x30 or 3-H-SB-4H-3x30	3-S-LL60 60cc Luer Lock Syringes	3-CTT-1H2-3x15 1hole, 3mm x 15cm



24H small fat parcels  
4H large fat parcels



1H2

1-800-525-7752

Cosmetic Surgery Suppliers Inc.

[www.cosmeticsurgerysuppliers.com](http://www.cosmeticsurgerysuppliers.com)

# ◆ Fat Transfer/Grafting ◆

*Simplify and Eliminate Many Timely and Trauma-Inducing Steps  
Involved in Current Autologous Fat Transfer Procedures*




## **Tissue-Trans® MEGA 1500**

Simple, sterile non-traumatic harvest of lipocytes.  
For all volume procedures:

- Eyelids
- Lips
- Face
- Hands
- Body Contouring
- Buttocks
- Breast Reconstruction
- Breast Augmentation
- Correction of Lipoplasty Deformities

The Tissue-Trans® does not change the normal flow of the physician's technique. Standard tumescent solution is infused. The Tissue-Trans® MEGA 1500 is placed between the harvest cannula and the waste canister. Fat is collected within the sterile Tissue-Trans® without exposure to ambient air. When harvest is complete, the MEGA 1500 is rotated and fat decanted. The waste (blood, epinephrine, Xylocaine, and plasma) is drained. Sterile fat is then drawn out with a luer lock syringe and is ready to inject to the recipient site.



Harvest up to 1500cc into sterile unit.  
Rotate unit to decant.  
Attach transfer tube and drain waste.  
Attach Luer lock syringes to draw out fat.  
Re-inject.

Instructional videos and white papers  
are available at:  
**[www.shippertmedical.com](http://www.shippertmedical.com)**



# Thermoplas Nasal Splint

## Nasal Splints



splints - 10 per package  
Actual Size

## Nasal Splint Kit



nasal splint kit - 5 kits per package

- Thermoplas is the first nasal splinting material with 100% memory, lets you undo what you have done & do it again perfectly.
- Thermoplas activates in ordinary hot tap water & needs no additional adhesives, no mess or fuss.
- Thermoplas is reasonably priced.
- Thermoplas is the first nasal splint that advances ease & patient comfort. There is no splint like it.

## Direction for use:

Place into hot water of at least 135°F for 15-20 seconds. Thermoplas will become transparent & moldable. Hardening will occur within approximately 60 seconds. To remold, simply place into hot water & it will return to original shape. If the Thermoplas splint does not become transparent & moldable, your hot water heater is not set at a high enough setting. You may also heat the Thermoplas splint by placing it in a container of water & using a microwave oven.

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)



# TLS® Drainage Systems

Let the **T**iny **L**ittle **S**ucker Do BIG Things for You!

## TLS Surgical Drainage System

For wound care of the extremities, the compact design of the TLS Surgical Drainage System is easily incorporated into the wound dressing and provides convenient and effective draining. This system utilizes standard 15cc evacuated blood collection tubes as reservoirs.

Packaged STERILE in a convenient tray, these kits are provided complete with radiopaque indicated, silicone catheter; trocar; protective housing; connecting adapter; and two 15cc evacuated tubes.

### CAT# DESCRIPTION

Kit contains: 1 Catheter, 1 Trocar, and 2 Evacuated Tubes

6640 7 FR, STERILE

6650 10 FR, STERILE

6660 4mm Flat with Side Holes, STERILE


The QUINTUBE® Monitor Pack, consisting of five additional evacuated tubes for purposes of monitoring drainage volume and color over time are available as replacement reservoirs for the TLS Surgical Drainage System.

### CAT# DESCRIPTION

6630 QUINTUBE Monitor Pack, 5 Tubes/Pack, 24 Packs/Case



## TLS Surgical Drainage System



## TLS Bulb Drainage System

## TLS Bulb Drainage System

For procedures requiring a larger reservoir, the TLS Bulb Drainage System provides an excellent alternative. These drainage systems provide gentle suction and prohibit reverse flow by relying on a 100cc silicone reservoir with an anti-reflux valve for constant, steady evacuation of postoperative fluids.

Packaged STERILE, these kits are provided complete with radiopaque indicated silicone catheter, trocar and reservoir.

### CAT# DESCRIPTION

Kit Contains: 1 Catheter, 1 Trocar, and 1 Reservoir

6646 15 FR, STERILE

6647 4mm Flat with Side Holes, STERILE

6648 10 FR, STERILE

6649 4mm Flat with Center Holes, STERILE

### Individual Components

6644 Trocar, 15 FR x 5"

6653 Catheter, 15 FR

6654 Silicone Bulb Reservoir, 100cc

## Ideal Kits for Hospital and Outpatient Surgical Centers

The TLS Surgical Drainage System and the TLS Bulb Drainage System are sold in sterile, convenient, complete kits consisting of a reservoir, catheter, and trocar. Both systems incorporate the small diameter silicone catheter and low profile trocars. Both systems provide positive suction drainage of fluid and tissue debris without reliance on gravity for a reduction in post operative complications and a decrease in patient healing time. The two systems permit accurate measurement and analysis of drainage fluids for optimal monitoring of both the wound and the patient.

# SQUEEZE-MARK® Surgical Markers

*Make Your Mark with SQUEEZE-MARK*




## **SQUEEZE-MARK Surgical Markers**

SQUEEZE-MARK Surgical Markers have a unique, patented, squeezable barrel that provides positive ink flow and a continuous marking action that allows marking over wet or prepped surfaces. The anti-rolling, ergonomic design combined with the squeezable barrel, makes SQUEEZE-MARK a popular choice among healthcare professionals.

The compact design of the SQUEEZE-MARK requires less inventory shelf space and results in less waste volume while providing the capabilities of a full length marker.

The SQUEEZE-MARK Surgical Marker is available in highly visible, scrub-resistant, gentian violet ink in both regular and micro tips. All SQUEEZE-MARK Surgical Markers are provided STERILE, packaged individually with a separate flexible ruler or adhesive ruler/labels.


CAT#	DESCRIPTION
9550	Regular Tip with Ruler, STERILE
9555	Regular Tip with Ruler/Labels, STERILE
9560	Micro Tip with Ruler, STERILE
9565	Micro Tip with Ruler/Labels, STERILE





# Endo-brow Kit

## Endoscope & Endoscope Sheath



Scope- 4mm x 15cm, 30 degrees, high magnification endoscope, autoclavable

Sheath- Visualization /irrigation sheath, open viewer spatula tip

## Rotational Cautery Instruments



micro metzenbaum scissor, curved



micro metzenbaum scissor, straight



micro straight jaw forcep



micro dissecting forcep, full curve tip



micro hook scissor

## Dissectors/Elevators



12mm broad dissector, chisel tip, curved

Periosteal elevator, 6mm oval, slight curved

Periosteal dissector, chisel tip, 7mm, slight curved

Periosteal elevator, spade end, 13mm, straight

Periosteal dissector, 7mm concave wedge, curved

Periosteal elevator 10mm, curved

3mm rounded micro dissector, slight curved

Nerve hook, tapered, curved

10mm round dissector, downturn tip, slight S-curve

Upturn tip, delicate dissector, slight S-curve

**A fully comprehensive endoplastic instrument set provided in a briefcase carrying case.**

**Set includes: A 30 degree endoscope and a fenestrated sheath; this shields the tissue while allowing full visualization and retraction. Includes built-in suction and irrigation port.**

**Five cautery handled, fully rotational instruments to adapt to left and right use without requiring separate instruments.**

**Ten of the most popular dissectors/elevators allowing for diverse choice and flexibility.**

**Significantly more cost effective than individual item purchases. This kit provides more than you need for far less than you would imagine.**

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# T.U.B.A. Set

---



## Set Includes:

- 1 - Trocar Dissector Handle
- 2 - Bullet Dissector Tips (1 small, 1 large)
- 2 - Trocar Tubes (1 small, 1 large)
- 1 - Hockey Stick Dissector

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.  
7758 Hampton Place  
Loganville, GA 30052  
800.525.7752  
770.934.0235  
[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)



# Hand Engine

---



35,000 RPM motor

Autoclavable hand piece and cord

Rotation easily switched from forward to reverse.

Self locking head - No more twisting and turning to remove. Just pull apart.

Use either manual speed control knob or variable speed foot switch.

Many accessories available

---

Purchase from:




Cosmetic Surgery Suppliers, Inc.  
7758 Hampton Place  
Loganville, GA 30052  
800.525.7752  
770.934.0235  
**cosmeticsurgerysuppliers.com**

# Diamond Fraises

## Diamond Fraises - for dermabrasion & cutaneous prophylaxis

The largest selection of sizes and grades of abrasion available anywhere

Actual Sizes



17x10 PEAR 17x8 17x15 15x10 DOME 17x6 BULLET 17x4 FRUSTUM 17x2.5 FLAME 13x2.5

Regular - Coarse - Extra Coarse

17x14mm.	\$65	\$124	\$145
17x10mm.	52	110	131
17x8mm.	44	94	111
15x10mm.	52	94	111
17x6mm.	37	79	93
17x4mm.	29	64	74
17x2.5mm.	24	47	57
13x2.5mm.	21	44	
10x2.5mm.	21	42	
CONE	23	30	
PEAR	23	30	37
DOME	91	123	160
BULLET	40	53	69
FRUSTUM	68	108	140
FLAME	32	43	52
BUTTON	34	46	59
MINI BULLET	26	36	53



**Top Row:** 10x2.5mm, Cone


**Bottom Row:** Button, Mini-Bullet

- one piece construction for 100% centricity
- diamond chips grade for grade are the largest and provide for maximum abrasion
- all edges have been rounded, avoiding dangerous sharp corners

## Stainless Steel Wire Planing Brushes

All fraises and brushes have 3/32" diameter shaft which is the international standard

Actual Sizes



Small - 17x3mm. (regular, soft, course)	\$69
Medium - 17x6mm.	72
Large - 17.9mm.	79
Extra Large - 17x12mm.	82

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Cosmetic Tattoo Kit

---



The REVOLUTION™ Cosmetic Kit consists of the following items:

- 1 Revolution Instrument with 3-needle stainless steel liner bar and round tip tube
  - 1 footswitch
  - 1 connecting cord
  - 10<sup>1/2</sup>-ounce bottles of cosmetic pigment  
(You may pick from our selection of colors)
  - 1 stainless steel round tip liner tube
- 1 stainless steel 4-needle flat square tip shader tube
- 1 package of five 4-needle S/S shader bars
- 1 package of five 3-needle S/S needle bars
- 1 package of five 1-needle S/S liner bars
  - 1 #16 cap holder for color
  - 1 package #16 caps
  - 1 screwdriver
- 1 power unit 2-power supply (110 volt)  
(Also available in 220 volt)

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place


Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Tattoo Pigment Chart



Purchase from:



Cosmetic Surgery Suppliers, Inc.  
 7758 Hampton Place  
 Loganville, GA 30052  
 800.525.7752  
 770.934.0235  
[cosmeticsurgesuppliers.com](http://cosmeticsurgesuppliers.com)



# Cosmetic Tattoo Price List

CATALOG	DESCRIPTION	PRICE
Sizes 1-8	Needle Tubes 1, 3, 4, 5, 6, 8	\$26.50
2000A	One needle liner bars - 5/pkg	\$24.50
2000B	Three needle liner bars - 5/pkg	\$24.50
2000F	Four needle liner bars - 5/pkg	\$27.50
2000D	Five needle liner bars - 5/pkg	\$26.50
2000G	Six needle shader bars - 5/pkg	\$33.50
2000H	Eight needle shader bars - 5/pkg	\$33.50
5023	#9 Cap Holder	\$4.00
5024	#16 Cap Holder	\$6.50
5026/5027	#9 Mixing Caps	\$4.00 - 18.00
5028/5029	#16 Mixing Caps	\$4.00 - 18.00

## COSMETIC COLORS

8001 Black	8019 Dusk #2
8002 Gray	8020 Sandstone #1
8003 Brown	8021 Sandstone #2
8004 Lt.Brown #1	8022 Flesh #1
8005 Lt.Brown #2	8023 Flesh #2
8006 Dark Brown	8024 Desert Brown #1
8007 White	8025 Desert Brown #2
8008 Cherri Berri #1	8026 Salmon
8009 Cherri Berri #2	8027 Tan #1
8010 Lemon Up #1	8028 Tan #2
8011 Lemon Up #2	8029 Sahara Tan #1
8012 Banana #1	8030 Sahara Tan #2
8013 Banana #2	8031 Indian Brown
8014 Peach #1	8032 Gobi Brown
8015 Peach #2	8033 Safari Tan #1
8016 Fire Red	8034 Safari Tan #2
8018 Dusk #1	8035 Black #1

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052


800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Skin Care Products

## Liposome CFT



4 oz. jar

Cosmetic Surgery Suppliers, Inc. is pleased to introduce a revolutionary product. Liposome CFT is a scientific breakthrough in controlling unsightly cellulite. Aminophylline, the key ingredient used in most slimming creams today, breaks down the fat so that it may be transported out of the body. Unlike other creams that have a 2% concentration, Liposome CFT has a 2.5%, the maximum concentration that will not cause skin irritation. What makes Liposome CFT different from the others is the added benefit of three other ingredients which increase the effectiveness & efficiency of the Aminophylline.

- L-Carnitine is a natural by-product of choline, which helps to prevent fatty deposits.
- Co-enzyme A is an important catalyst in penetrating the skin & fatty tissue.
- Aloe, being the primary ingredient in Liposome CFT, (as opposed to water as in the other similar products), plus Glycerin and other natural extracts makes this slimming cream soothing to the skin, as well as moisturizing.

## Dermal-K

A Vitamin K Cream that promotes healing and reduces bruising this cream utilizes a unique all natural dermaceytical delivery system that is now available from Cosmetic Surgery Suppliers, Inc.

Used to treat spider veins, promote normal healing, reduce bruising, relieve pain and inflammation in wounds. Dermal K, also may be helpful in treating burns, surface ulcers, cuts, abrasions, and skin irritation. The all natural ingredients have a synergistic effect and it is this combination which stimulates epithelial growth and promotes healthy granulation and clarifying of injured skin.

The use of Vitamin K injections to treat bleeding problems is well documented and now the same benefits are available in the form of an all natural topical cream. Used pre-operatively, twice a day, for seven days, Dermal-K has been shown to effectively reduce post-op bruising.



2 oz. jar

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# The Real Image Mirror™

---



Have you seen yourself the way others see you?

---

## The Real Image Mirror™

- A positive real image of the patient - reversed from standard mirrors.
- Allows patients to see themselves as they are seen by others.
- Allows patient to see their facial asymmetries.
- Prevents an "unhappy" patient post-operatively from thinking YOU caused their naturally existing asymmetry.
- Hand-held and portable.
- "First surface" mirror eliminates distortion of the reflection.

---

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

[cosmeticsurgerysuppliers.com](http://cosmeticsurgerysuppliers.com)

# Swiss Therapy™

## Masks and Pads

Invotec's Swiss Therapy™ Masks and Pads are used as a lightweight and flexible cold compress during and after cosmetic procedures of the face. Synthesized in Switzerland of a patented Hydrogel material, Swiss Therapy is similar to natural tissue, making the product comfortable and highly compliant.

### **Blepharoplasty • Rhinoplasty • Laser Skin Resurfacing • Dermabrasion**

- *Moist, but not wet*
- *Lightweight & very comfortable*
- *Cools by contact & evaporation*
- *Gently pliable and conforming*
- *Non-adherent*
- *Highly patient compliant*
- *Convenient, easy to use*
- *Sterile - 2 year shelf life*



#### **Swiss Therapy™ Masks**

Swiss Therapy™ Masks, used after cosmetic procedures, are a cooling PostOP compress for eye lid surgery such as Blepharoplasty. While comforting to the patient, it reduces redness and swelling, increasing the patient's satisfaction.

*Gentle      Comfortable      Cooling*

#### **LASER & IPL Skin Therapy with Swiss Therapy™**

Use Swiss Therapy™ pads during skin resurfacing with Laser and IPL<sup>1</sup> to ensure even cooling of a large area and prevent erythemas formation. Optimize the treatment period and therapeutic outcome with Swiss Therapy from Invotec.

*Increased epidermal transparency  
Improved visibility of skin structure  
Better penetration of the site  
Reduces scatter of laser light  
Less than 1% power loss  
Soothing to skin*



Catalog No.	Description		Qty/Box
10-08803	Swiss Therapy™ Eye Mask	Sterile	3
10-08805	Swiss Therapy™ Sheet 100 x 100 x 3.3mm	Sterile	5

#### **To order, please contact:**



**COSMETIC SURGERY SUPPLIERS**  
**7758 HAMPTON PLACE**  
**LOGANVILLE, GA 30052**  
 PH (770) 934-0235 • (800) 525-7752  
 FAX (770) 939-5755 • (800) 574-3563

*Cooling  
Contact  
Transparency*

<sup>1</sup>Swiss Therapy Pads can be used with all laser systems (Q-switched, pulsed) and arc lamps in a wavelength range of 400-1100nm without affecting transmission of the laser energy. Swiss Therapy should not be used with wavelengths less than 400nm or more than 1100nm or in laser therapies using lasers with integrated, liquid coolant based cooling systems.