

Tumescent Infiltration

The Fluid Transfer Kit

Kit contains:

- One sterile transfer system with dual-check valve and bag spike
- One sterile 60cc syringe
- Priced lower than current systems

For use where ever closed fluid transfer is desired. Large 60cc syringe allows rapid movement of fluids with fewer repetitions. Dual-check valve means no disconnect/reconnect problems.



Blood Pressure Cuff

- Pressure infuser
- Fast fluid level inspection
- Durable and easy to clean
- Precise easily read pressure gauge



Variable Speed Pump

- Flow rate controlled and adjustable from 35 - 350 ml/min
- Easy set up. Disposable tubing set
- Continuous infusion
- Pneumatic footswitch included
- Directional switch - tubing can be installed in either direction



Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

cosmeticsurgesuppliers.com

Tumescent Instruments and Tubing



Infusion Handle



Infusion Handle with On/Off Control



Tumescent Needle

Luer Lock Connector With Multi-Directional Holes For Rapid Infusion of
Tumescent Fluids

Standard Sizes

2mm x 15cm

2mm x 23cm

2mm x 30cm

2.5mm x 15cm

2.5mm x 23cm

2.5mm x 30cm

3mm x 15cm

3mm x 23cm

3mm x 30cm



AFT Tubing

(For use with electric pumps)



T.I. Tubing

(For use with pressure bag)

Purchase from:



Cosmetic Surgery Suppliers, Inc.

7758 Hampton Place

Loganville, GA 30052

800.525.7752

770.934.0235

cosmeticsurgerysuppliers.com